

Recombinant Human Beta-defensin-2 (BD-2)

CATALOG #:	4372-20	20 µg
	4372-1000	1 mg
LOT #:	_____	
SOURCE:	<i>E. coli</i>	
PURITY:	>98% by SDS-PAGE and HPLC analyses Endotoxin level is <0.1 ng per µg (1EU/µg).	
MOL. WEIGHT:	4.3 kDa	
FORM:	Sterile filtered & Lyophilized with no additives	

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in ddH₂O to a concentration of 0.1-1.0 mg/ml. The solution can then be diluted into other aqueous buffers or stored at 4°C for 1 week or -20°C for future use.

STORAGE CONDITIONS:

The lyophilized protein is best stored desiccated below 0°C. Reconstituted human BD-2 should be stored at working aliquots at -20°C.

BACKGROUND DESCRIPTION:

Recombinant human beta-defensin-2 (BD-2) is an antimicrobial peptide which belongs to the distinct family of beta-defensins. Recombinant human BD-2 is a 4.3 kDa protein containing 41 amino acid residues.

BIOLOGICAL ACTIVITY:

Determined by its ability to chemoattract immature human dendritic cells using a concentration of 10-100 ng/ml.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used on humans.