

Recombinant Human Visfatin

CATALOG #:	4907-10	10 µg
	4907-50	50 µg
	4907-1000	1 mg
LOT #:		
SOURCE:	Produced in <i>E. coli</i> . Full length human visfatin (a.a. 1-491) is fused at the N-terminus to a His-Tag.	
MOL. WT.	55 kDa	
PURITY:	>95% by SDS-PAGE	
ENDOTOXIN:	≤0.1 EU/µg protein (LAL test)	
FORM:	Lyophilized with no additives.	
RECONSTITUTION:	Reconstitute in PBS to the concentration of 0.1-1 mg/ml.	
SHIPPING:	Ship on blue ice	
STORAGE:	-20°C	
DESCRIPTION:	Pre-B-Cell Colony Enhancing Factor 1	
HANDLING:	For maximum product recovery, centrifuge the vial prior to opening. After reconstitution in PBS, aliquot and store at -20°C. Stable for up to 3 months. Avoid freeze/thaw cycle.	
BIOLOGICAL ACTIVITY:	Mimics insulin effects in differentiating 3T3-L1 cells.	

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

FOR RESEARCH USE ONLY! Not to be used on humans.