

**Catalog Number:** BDH83634.400

**Catalog Description:** N,N-Dimethylformamide

**Grade:** HiPerSolv CHROMANORM for HPLC

**Date Current Specification Set:** 2-16-16

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

### Characteristics

### Specifications

Assay (by GC-FID)

Min 99.90 %

Water

(by KF coulometry)

Max 0.02 %

Nonvolatile Residue

Max 0.0005 %

Acidity

Max 0.0005 meq/g

Alkalinity

Max 0.0002 meq/g

Transmittance at

270 nm

Min 10.0 %

275 nm

Min 50.0 %

290 nm

Min 80.0 %

300 nm


Min 90.0 %

330 nm

Min 98.0 %

Signed on behalf of VWR:

Quality Control Manager



**Additional Information (optional):**