

Product Specifications

	VWR Part Number:	IR2717
	ImmunoReagents, Inc. Part Number:	MAB-100-IL6
Product:	MONOCLONAL ANTI-IL-6	
Concentration:	1 mg/mL solution	
Amount:	100 µg	
Conjugate:	Unconjugated	
Form:	1 mg/mL solution	
Purification:	Protein A	
Purity:	> 95% based on SDS-PAGE	
Host:	Mouse	
Immunogen:	Interleukin-6	
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with no stabilizers or preservatives	
Storage:	Store at 2-8°C for use within 1-3 weeks or -20°C for long term storage, avoid repeat freeze thaws.	
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.	
Specificity:	Interleukin-6	
Cross Reactivity:	Interleukin 6 (IL-6) is an important pro-inflammatory and anti-inflammatory cytokine expressed by T cells, macrophages and muscle cells. IL-6 signals through a receptor complex containing two receptors, IL-6R α and gp130. IL-6 has an important function in	
Country of Origin:	US Origin	
Applications:	FC	
Disclaimer:	For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.	

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.