

Product Specifications

	VWR Part Number:	IR2723
	ImmunoReagents, Inc. Part Number:	MAB-200-P28BIO
Product:	MONOCLONAL ANTI-IL-27 P28 SUBUNIT	
Concentration:	0.5 mg/mL solution	
Amount:	100 µg	
Conjugate:	Biotin (NHS-LC)	
Form:	0.5 mg/mL solution	
Purification:	Protein A	
Purity:	> 95% based on SDS-PAGE	
Host:	Mouse	
Immunogen:	Interleukin-27 p28 Subunit	
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with no stabilizers or preservatives	
Storage:	Store at 2-8°C for use within 1-3 weeks or -20°C for long term storage, avoid repeat freeze thaws.	
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.	
Specificity:	Interleukin-27 p28 Subunit	
Cross Reactivity:	The p28 subunit of IL-27 (IL-27/p28), also known as Interleukin-30 (IL-30), is a member of the IL-12 family of cytokines. When combined with EB13 (Epstein-Barr virus induced gene 3), p28 forms a heterodimer known as IL-27. p28 is thought to be a pro-inf	
Country of Origin:	US Origin	
Applications:	FC	
Disclaimer:	For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.	

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.