
Product Specifications

	VWR Part Number:	IR2765
	ImmunoReagents, Inc. Part Number:	PAR-200-MIP3A
Product:	RABBIT ANTI-MIP-3A	
Concentration:	N/A	
Amount:	100 µg	
Conjugate:	Unconjugated	
Form:	Lyophilized	
Purification:	Protein A	
Purity:	> 90% based on SDS-PAGE	
Host:	Rabbit	
Immunogen:	Mouse MIP-3?	
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 with no stabilizers or preservatives	
Storage:	Store at 2-8°C after restoration with 0.1mL deionized water or -20°C for long term storage, avoid repeat freeze thaws.	
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.	
Specificity:	Mouse MIP-3?	
Cross Reactivity:	N/A	
Country of Origin:	US Origin	
Applications:	Western blot analysis of recombinant mouse RANK-L (200-04) using rabbit anti- mouse RANK-L polyclonal antibody and anti-rabbit secondary antibodies. MW of mouse RANK-L is ~20 kDa.	
Disclaimer:	For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not for human or animal consumption. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license under any patent of ImmunoReagents, Inc.	

Product may not be resold or modified for resale without prior written approval of ImmunoReagents, Inc.