

## Product Specifications

	<b>VWR Part Number:</b>	<b>DASHRP-015</b>
	<b>ImmunoReagents, Inc. Part Number:</b>	<b>DASHRP-015</b>
<b>Product:</b>	Ultra-Poly HRP Anti-Sheep IgG (H&L)	
<b>Concentration:</b>	N/A	
<b>Amount:</b>	15ml	
<b>Conjugate:</b>	Horseradish Peroxidase	
<b>Form:</b>	Liquid	
<b>Purification:</b>	Liquid	
<b>Purity:</b>	Affinity purified antibody is > 95% based on SDS-PAGE	
<b>Host:</b>	N/A	
<b>Immunogen:</b>	Sheep IgG	
<b>Buffer:</b>	PBS & 0.1% proclin 150	
<b>Storage:</b>	2-8 °C	
<b>Shelf Life:</b>	1 year from date of receipt. Prepare working dilution prior to use and then discard.	
<b>Specificity:</b>	Based on IEP, this immunopolymer antibody reacts with heavy (?) chains on sheep IgG, light chains on all sheep immunoglobulins	
<b>Cross Reactivity:</b>	Based on IEP, no cross-reactivity is observed to non-immunoglobulin sheep serum proteins, IgG from human, mouse or rabbit	
<b>Country of Origin:</b>	US Origin	
<b>Applications:</b>	IHC	
<b>Disclaimer:</b>	N/A	