

Product Specifications

	VWR Part Number:	DASHRP-050
	ImmunoReagents, Inc. Part Number:	DASHRP-050
Product:	Ultra-Poly HRP Anti-Sheep IgG (H&L)	
Concentration:	N/A	
Amount:	50ml	
Conjugate:	Horseradish Peroxidase	
Form:	Liquid	
Purification:	Liquid	
Purity:	Affinity purified antibody is > 95% based on SDS-PAGE	
Host:	N/A	
Immunogen:	Sheep IgG	
Buffer:	PBS & 0.1% proclin 150	
Storage:	2-8 °C	
Shelf Life:	1 year from date of receipt. Prepare working dilution prior to use and then discard.	
Specificity:	Based on IEP, this immunopolymer antibody reacts with heavy (?) chains on sheep IgG, light chains on all sheep immunoglobulins	
Cross Reactivity:	Based on IEP, no cross-reactivity is observed to non-immunoglobulin sheep serum proteins, IgG from human, mouse or rabbit	
Country of Origin:	US Origin	
Applications:	IHC	
Disclaimer:	N/A	