

Revision date: 06-05-2014

# SAFETY DATA SHEET

# 1. Identification

**Product identifier: POTASSIUM HYDROXIDE SOLUTION** 

Other means of identification

Product No.: 3170, 3149, 3147, 3144, 3143, 6671, 3148, 73317, 37828, 23600

Recommended use and restriction on use

Recommended use: Not available. Restrictions on use: Not known.

# Manufacturer/Importer/Supplier/Distributor information

Manufacturer

Company Name: Avantor Performance Materials, Inc. Address: 3477 Corporate Parkway, Suite 200

Center Valley, PA 18034

Telephone:

Customer Service: 855-282-6867

Fax:

Contact Person: Environmental Health & Safety e-mail: info@avantormaterials.com

**Emergency telephone number:** 

24 Hour Emergency: 908-859-2151

Chemtrec: 800-424-9300

# 2. Hazard(s) identification

#### Hazard classification

## Physical hazards

Corrosive to metals Category 1

**Health hazards** 

Skin corrosion/irritation Category 1
Specific target organ toxicity - single Category 3

exposure

# Unknown toxicity

Acute toxicity, oral 0 %
Acute toxicity, dermal 43 %
Acute toxicity, inhalation, vapor 100 %
Acute toxicity, inhalation, dust or mist 100 %

**Unknown toxicity** 

Chronic hazards to the aquatic 43 %

environment

## Label elements

# Hazard symbol:



Revision date: 06-05-2014



Signal word: Danger

**Hazard statement:** May be corrosive to metals.

Harmful if swallowed.

Causes severe skin burns and eye damage.

May cause respiratory irritation. May cause drowsiness or dizziness.

Harmful to aquatic life.

# **Precautionary statement**

**Prevention:** Do not breathe dust/fume/gas/mist/vapors/spray. Keep only in original

container. Use only outdoors or in a well-ventilated area. Wash thoroughly

after handling. Wear protective gloves/protective clothing/eye

protection/face protection. Do not eat, drink or smoke when using this product. Avoid breathing dust/fume/gas/mist/vapors/spray. Avoid release to

the environment.

**Response:** Absorb spillage to prevent material damage. IF SWALLOWED: Rinse

mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON

CENTER or doctor/physician.

Storage: Store locked up. Store in a well-ventilated place. Keep container tightly

closed. Store in corrosive resistant container with a resistant inner liner.

**Disposal:** Dispose of contents/container to an appropriate treatment and disposal

facility in accordance with applicable laws and regulations, and product

characteristics at time of disposal.

Other hazards which do not result in GHS classification:

None.

# 3. Composition/information on ingredients

# Mixtures

Chemical identity	Common name and synonyms	CAS number	Content in percent (%)*
POTASSIUM HYDROXIDE		1310-58-3	40 - 60%

<sup>\*</sup> All concentrations are percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.

## 4. First-aid measures

General information: Get medical advice/attention if you feel unwell. Show this safety data sheet

to the doctor in attendance.



Revision date: 06-05-2014

**Ingestion:** Call a physician or poison control center immediately. Rinse mouth. Do

NOT induce vomiting. If vomiting occurs, keep head low so that stomach

content doesn't get into the lungs.

**Inhalation:** Move to fresh air. Call a physician or poison control center immediately. If

breathing is difficult, give oxygen. Apply artificial respiration if victim is not

breathing

**Skin contact:** Immediately flush with plenty of water for at least 15 minutes while

removing contaminated clothing and shoes. Call a physician or poison control center immediately. Wash contaminated clothing before reuse.

Destroy or thoroughly clean contaminated shoes.

**Eye contact:** Immediately flush with plenty of water for at least 15 minutes. If easy to do,

remove contact lenses. Call a physician or poison control center

immediately. In case of irritation from airborne exposure, move to fresh air.

Most important symptoms/effects, acute and delayed

**Symptoms:** Causes severe skin and eye burns. Causes digestive tract burns. Mist or

vapor extremely irritating to eyes and respiratory tract.

Indication of immediate medical attention and special treatment needed

**Treat symptomatically.** Symptoms may be delayed.

5. Fire-fighting measures

General fire hazards: No data available.

Suitable (and unsuitable) extinguishing media

Suitable extinguishing

media:

Water spray, foam, dry powder or carbon dioxide.

Unsuitable extinguishing

media:

None known.

Specific hazards arising from

the chemical:

During fire, gases hazardous to health may be formed. Wear appropriate protective gear if spilled during fire fighting. Product is highly caustic.

Special protective equipment and precautions for firefighters

Special fire fighting

procedures:

Move containers from fire area if you can do so without risk. Use water spray to keep fire-exposed containers cool. Cool containers exposed to

flames with water until well after the fire is out.

Special protective equipment

for fire-fighters:

Firefighters must use standard protective equipment including flame retardant coat, helmet with face shield, gloves, rubber boots, and in

enclosed spaces, SCBA.

## 6. Accidental release measures

Personal precautions, protective equipment and emergency procedures: Keep unauthorized personnel away. Keep upwind. Ventilate closed spaces before entering them. Use personal protective equipment. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. See Section 8 of the MSDS for Personal Protective Equipment.



Revision date: 06-05-2014

Methods and material for containment and cleaning

up:

Neutralize spill area and washings with dilute acetic acid. Absorb spill with vermiculite or other inert material, then place in a container for chemical waste. Dike far ahead of larger spill for later recovery and disposal.

Notification Procedures: Dike for later disposal. Prevent entry into waterways, sewer, basements or

confined areas. Stop the flow of material, if this is without risk. Inform

authorities if large amounts are involved.

**Environmental precautions:** Do not contaminate water sources or sewer. Prevent further leakage or

spillage if safe to do so. Avoid discharge into drains, water courses or onto

the ground.

# 7. Handling and storage

**Precautions for safe handling:** Wear protective gloves/protective clothing/eye protection/face protection.

Do not get in eyes, on skin, on clothing. Do not taste or swallow. Do not breathe dust/fume/gas/mist/vapors/spray. Wash hands thoroughly after handling. Do not eat, drink or smoke when using the product. See Section 8 of the MSDS for Personal Protective Equipment. Use caution when adding

this material to water.

Conditions for safe storage,

including any incompatibilities:

Keep away from food, drink and animal feeding stuffs. Keep container tightly closed. Do not store in metal containers. Store in a cool and well-

ventilated place. Store in a dry place.

# 8. Exposure controls/personal protection

## **Control parameters**

Occupational exposure limits

Chemical identity	Туре	Exposure Limit values	Source
POTASSIUM HYDROXIDE	Ceiling	2 mg/m3	US. ACGIH Threshold Limit Values (2011)
	REL	2 mg/m3	US. NIOSH: Pocket Guide to Chemical Hazards (2010)
	Ceiling	2 mg/m3	US. OSHA Table Z-1-A (29 CFR 1910.1000) (1989)

# Appropriate engineering controls

No data available.

#### Individual protection measures, such as personal protective equipment

**General information:** Good general ventilation (typically 10 air changes per hour) should be used.

Ventilation rates should be matched to conditions. If applicable, use

process enclosures, local exhaust ventilation, or other engineering controls to maintain airborne levels below recommended exposure limits. If exposure limits have not been established, maintain airborne levels to an acceptable level. An eye wash and safety shower must be available in the

immediate work area.

**Eye/face protection:** Wear safety glasses with side shields (or goggles) and a face shield. Wear

a full-face respirator, if needed.

Skin protection

**Hand protection:** Chemical resistant gloves

**Other:** Wear suitable protective clothing.



Revision date: 06-05-2014

**Respiratory protection:** If engineering controls do not maintain airborne concentrations below

recommended exposure limits (where applicable) or to an acceptable level (in countries where exposure limits have not been established), an approved respirator must be worn. Air-purifying respirator with an appropriate, government approved (where applicable), air-purifying filter,

cartridge or canister. Contact health and safety professional or

manufacturer for specific information.

**Hygiene measures:** Provide eyewash station and safety shower. Always observe good personal

hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking. Routinely wash work clothing and

protective equipment to remove contaminants. Do not get in eyes. Do not

get this material in contact with skin.

# 9. Physical and chemical properties

**Appearance** 

Physical state: Liquid
Form: Liquid
Color: Colorless
Odor: Odorless

Odor threshold: No data available.

pH: 13.5 0.1 N Aqueous solution

Melting point/freezing point: -29 °C Initial boiling point and boiling range: 132 °C

Flash Point: Not applicable Evaporation rate: As water

Flammability (solid, gas): No data available.

Upper/lower limit on flammability or explosive limits

Flammability limit - upper (%):

Flammability limit - lower (%):

Explosive limit - upper (%):

No data available.

No data available.

No data available.

No data available.

Vapor pressure: 0.27 kPa

**Vapor density:**No data available. **Relative density:**1.45 (20 °C)

Solubility(ies)

Solubility in water: Soluble

Solubility (other):

Partition coefficient (n-octanol/water):

Auto-ignition temperature:

Decomposition temperature:

Viscosity:

No data available.

No data available.

No data available.

No data available.

# 10. Stability and reactivity

**Reactivity:** Reacts violently with strong acids. May react with strong oxidizers.

**Chemical stability:** Material is stable under normal conditions.

Possibility of hazardous

reactions:

Hazardous polymerization does not occur.

**Conditions to avoid:** Excessive heat. Contact with incompatible materials.



Revision date: 06-05-2014

**Incompatible materials:** Acids. Strong oxidizing agents. Contact with metals may evolve flammable

hydrogen gas.

**Hazardous decomposition** 

products:

Oxides of potassium.

# 11. Toxicological information

# Information on likely routes of exposure

**Ingestion:** May cause burns of the gastrointestinal tract if swallowed.

**Inhalation:** May cause damage to mucous membranes in nose, throat, lungs and

bronchial system.

**Skin contact:** Causes severe skin burns.

**Eye contact:** Causes serious eye damage.

# Information on toxicological effects

## Acute toxicity (list all possible routes of exposure)

Oral

**Product:** ATEmix (Rat): 580.85 mg/kg

**Dermal** 

**Product:** 

No data available.

Inhalation

**Product:** No data available.

Repeated dose toxicity

**Product:** No data available.

Skin corrosion/irritation

**Product:** Causes severe skin burns.

Serious eye damage/eye irritation

**Product:** Causes serious eye damage.

Respiratory or skin sensitization

**Product:** Not a skin sensitizer.

Carcinogenicity

**Product:** This substance has no evidence of carcinogenic properties.

# IARC Monographs on the Evaluation of Carcinogenic Risks to Humans:

No carcinogenic components identified

# **US. National Toxicology Program (NTP) Report on Carcinogens:**

No carcinogenic components identified

#### US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050):

No carcinogenic components identified



Revision date: 06-05-2014

# Germ cell mutagenicity

In vitro

**Product:** No mutagenic components identified

In vivo

**Product:** No mutagenic components identified

Reproductive toxicity

**Product:** No components toxic to reproduction

Specific target organ toxicity - single exposure

**Product:** Respiratory tract irritation.

Specific target organ toxicity - repeated exposure

**Product:** None known.

**Aspiration hazard** 

**Product:** Not classified

Other effects: None known.

# 12. Ecological information

# **Ecotoxicity:**

## Acute hazards to the aquatic environment:

Fish

**Product:** No data available.

Specified substance(s):

POTASSIUM LC 50 (Western mosquitofish (Gambusia affinis), 96 h): 80 mg/l Mortality

**HYDROXIDE** 

**Aquatic invertebrates** 

**Product:** No data available.

## Chronic hazards to the aquatic environment:

**Fish** 

**Product:** No data available.

Aquatic invertebrates

**Product:** No data available.

**Toxicity to Aquatic Plants** 

**Product:** No data available.

# Persistence and degradability

**Biodegradation** 

**Product:** Expected to be readily biodegradable.

**BOD/COD** ratio

**Product:** No data available.

# **Bioaccumulative potential**

**Bioconcentration factor (BCF)** 

**Product:** No data available on bioaccumulation.

# Partition coefficient n-octanol / water (log Kow)

**Product:** No data available.



Revision date: 06-05-2014

**Mobility in soil:** The product is water soluble and may spread in water systems.

Other adverse effects: Expected to be harmful to aquatic organisms. The product may affect the

acidity (pH-factor) in water with risk of harmful effects to aquatic organisms.

# 13. Disposal considerations

**Disposal instructions:** Discharge, treatment, or disposal may be subject to national, state, or local

laws.

**Contaminated packaging:** Since emptied containers retain product residue, follow label warnings even

after container is emptied.

# 14. Transport information

DOT

UN number: UN 1814

UN proper shipping name: Potassium hydroxide, solution

Transport hazard class(es)

Class(es): 8
Label(s): 8
Packing group: II
Marine Pollutant: No

**IMDG** 

UN number: UN 1814

UN proper shipping name: POTASSIUM HYDROXIDE SOLUTION

Transport hazard class(es)

Class(es): 8 Label(s): 8

EmS No.: F-A, S-B

Packing group: II
Marine Pollutant: No

**IATA** 

UN number: UN 1814

Proper Shipping Name: Potassium hydroxide solution

Transport hazard class(es):

Class(es): 8
Label(s): 8

Marine Pollutant: No
Packing group: II

# 15. Regulatory information

# **US federal regulations**

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

None present or none present in regulated quantities.

CERCLA Hazardous Substance List (40 CFR 302.4):

POTASSIUM HYDROXIDE Reportable quantity: 1000 lbs.



Revision date: 06-05-2014

# Superfund amendments and reauthorization act of 1986 (SARA)

# **Hazard categories** Acute (Immediate) Chronic (Delayed) Fire Pressure Generating Reactive SARA 302 Extremely hazardous substance None present or none present in regulated quantities. SARA 304 Emergency release notification **Chemical identity** 1000 lbs. **POTASSIUM HYDROXIDE** SARA 311/312 Hazardous chemical **Threshold Planning Quantity** Chemical identity **POTASSIUM** 500 lbs **HYDROXIDE** SARA 313 (TRI reporting) None present or none present in regulated quantities.

## Clean Water Act Section 311 Hazardous Substances (40 CFR 117.3)

POTASSIUM HYDROXIDE Reportable quantity: 1000 lbs.

# Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130):

None present or none present in regulated quantities.

# **US state regulations**

# **US. California Proposition 65**

No ingredient regulated by CA Prop 65 present.

## US. New Jersey Worker and Community Right-to-Know Act

POTASSIUM HYDROXIDE Listed

#### **US. Massachusetts RTK - Substance List**

**POTASSIUM** Listed

**HYDROXIDE** 

# US. Pennsylvania RTK - Hazardous Substances

POTASSIUM HYDROXIDE Listed

## US. Rhode Island RTK

POTASSIUM HYDROXIDE Listed



Revision date: 06-05-2014

# **Inventory Status:**

Australia AICS: Canada DSL Inventory List: EINECS, ELINCS or NLP:

Japan (ENCS) List:

China Inv. Existing Chemical Substances: Korea Existing Chemicals Inv. (KECI):

Canada NDSL Inventory: Philippines PICCS: US TSCA Inventory:

New Zealand Inventory of Chemicals:


Japan ISHL Listing:

Japan Pharmacopoeia Listing:

On or in compliance with the inventory Not in compliance with the inventory. On or in compliance with the inventory Not in compliance with the inventory. On or in compliance with the inventory On or in compliance with the inventory On or in compliance with the inventory Not in compliance with the inventory. Not in compliance with the inventory.

# 16.Other information, including date of preparation or last revision

# **NFPA Hazard ID**



Hazard rating: 0 - Minimal; 1 - Slight; 2 - Moderate; 3 - Serious; 4 - Severe

COR: Corrosive

**Issue date:** 06-05-2014

**Revision date:** No data available.

Version #: 1.0

Further information: No data available.



Revision

Revision date: 06-05-2014

Version: 1.0

#### Disclaimer:

THE INFORMATION PRESENTED IN THIS MATERIAL SAFETY DATA SHEET (MSDS/SDS) WAS PREPARED BY TECHNICAL PERSONNEL BASED ON DATA THAT THEY BELIEVE IN THEIR GOOD FAITH JUDGMENT IS ACCURATE. HOWEVER, THE INFORMATION PROVIDED HEREIN IS PROVIDED "AS IS," AND AVANTOR PERFORMANCE MATERIALS MAKES AND GIVES NO REPRESENTATIONS OR WARRANTIES WHATSOEVER, AND EXPRESSLY DISCLAIMS ALL WARRANTIES REGARDING SUCH INFORMATION AND THE PRODUCT TO WHICH IT RELATES, WHETHER EXPRESS, IMPLIED, OR STATUTORY, INCLUDING WITHOUT LIMITATION, WARRANTIES OF ACCURACY, COMPLETENESS, MERCHANTABILITY, NON-INFRINGEMENT, PERFORMANCE, SAFETY, SUITABILITY, STABILITY, AND FITNESS FOR A PARTICULAR PURPOSE, AND ANY WARRANTIES ARISING FROM COURSE OF DEALING, COURSE OF PERFORMANCE. OR USAGE OF TRADE. THIS MSDS/SDS IS INTENDED ONLY AS A GUIDE TO THE APPROPRIATE PRECAUTIONARY HANDLING OF THE MATERIAL BY A PROPERLY TRAINED PERSON USING THIS PRODUCT, AND IS NOT INTENDED TO BE COMPREHENSIVE AS TO THE MANNER AND CONDITIONS OF USE, HANDLING, STORAGE, OR DISPOSAL OF THE PRODUCT. INDIVIDUALS RECEIVING THIS MSDS/SDS MUST ALWAYS EXERCISE THEIR OWN INDEPENDENT JUDGMENT IN DETERMINING THE APPROPRIATENESS OF SUCH ISSUES. ACCORDINGLY, AVANTOR PERFORMANCE MATERIALS ASSUMES NO LIABILITY WHATSOEVER FOR THE USE OF OR RELIANCE UPON THIS INFORMATION. NO SUGGESTIONS FOR USE ARE INTENDED AS, AND NOTHING HEREIN SHALL BE CONSTRUED AS, A RECOMMENDATION TO INFRINGE ANY EXISTING PATENTS OR TO VIOLATE ANY FEDERAL, STATE, LOCAL, OR FOREIGN LAWS. AVANTOR PERFORMANCE MATERIALS REMINDS YOU THAT IT IS YOUR LEGAL DUTY TO MAKE ALL INFORMATION IN THIS MSDS/SDS AVAILABLE TO YOUR EMPLOYEES.