

KIMAX® STARTER PACKS

• SAVE TIME • SAVE SPACE • SAVE MONEY WITH KIMAX® LAB STARTER PACKS

- Purchase only what you need
- Less storage space required
- The most popular sizes in one case saves time searching for the correct glassware



LAB STARTER PACKS



To Deliver Graduated Cylinder Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	±Tol. (mL)	Approx. Ht. (mm)	Case Qty.
89000-968	10	1 to 10	0.2	0.1	135	1
	25	3 to 25	0.5	0.3	140	1
	50	3 to 50	1	0.4	190	1
	100	5 to 100	1	0.6	255	1
	250	10 to 250	2	1.4	330	1



Low-Form Beaker Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Approx. OD x Ht. (mm)	Case Qty.
89000-924	50	10 to 40	10	42 x 53	1
	100	20 to 80	10	50 x 70	1
	250	25 to 200	10	68 x 88	1
	600	50 to 500	25	88 x 122	1
	1000	50 to 950	50	107 x 145	1



Heavy-Duty Beaker Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Approx. OD x Ht. (mm)	Case Qty.
89209-302	250	25 to 200	25	67 x 92	1
	400	25 to 325	50	77 x 114	1
	600	50 to 500	50	87 x 122	1
	1000	100 to 1000	50	107 x 152	1



Erlenmeyer Flask Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Widest Dia. (mm)	Approx. Ht. (mm)	Case Qty.
89001-342	50	20 to 50	10	50	78	1
	125	50 to 125	25	66	108	1
	250	50 to 225	25	81	130	1
	500	100 to 500	50	102	174	1
	1000	250 to 1000	50	128	213	1



Media Bottles Starter Pack

Catalog No.	Cap. (mL)	Grad. Interval (mL)	Subdiv. (mL)	Widest Dia. (mm)	Approx. Ht. (mm)	Case Qty.
10798-966	50	20 to 50	10	50	78	1
	125	50 to 125	25	66	108	1
	250	50 to 225	25	81	130	1
	500	100 to 500	50	102	174	1
	1000	250 to 1000	50	128	213	1



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.

1.800.932.5000 | vwr.com