



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Salmonella Antibody [6371]

CATALOG NUMBER: 35-643

Specifications

SPECIES REACTIVITY:

TESTED APPLICATIONS:

USER NOTE: Optimal dilutions for each application to be determined by the researcher.

SPECIFICITY: Specific for common core of all Salmonellae tested.

HOST SPECIES: Mouse

Properties

PURIFICATION: Salmonella antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.

PHYSICAL STATE: Liquid

BUFFER: Salmonella sp. antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

STORAGE CONDITIONS: Salmonella antibody can be stored at -20°C, stable for one year.

CLONALITY: Monoclonal

ISOTYPE: Mouse IgG2b

CONJUGATE: Unconjugated

Additional Info

ALTERNATE NAMES: Anti-Salmonella

Background

FOR RESEARCH USE ONLY

December 12, 2016