



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064



Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Trichomonas vaginalis Antibody [7111]

CATALOG NUMBER: 35-661

Specifications

USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Specific to the p65 adhesive antigen of <i>T. vaginalis</i> .
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Trichomonas vaginalis antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
PHYSICAL STATE:	Liquid
BUFFER:	Trichomonas vaginalis antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
STORAGE CONDITIONS:	Trichomonas vaginalis antibody can be stored at -20°C, stable for one year.
CLONALITY:	Monoclonal
ISOTYPE:	Mouse IgG1
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	Anti-Trichomonas vaginalis
-------------------------	----------------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016