



prosci-inc.com



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064

HIGH PERFORMANCE ANTIBODIES ... AND MORE

Toll Free: +1 (888) 513 9525  
Local: +1 (858) 513 2638  
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

## Verotoxins (Shiga-like) Antibody [7585]

CATALOG NUMBER: 35-675

### Specifications

<b>SPECIES REACTIVITY:</b>	Bacteria
<b>TESTED APPLICATIONS:</b>	ELISA
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>SPECIFICITY:</b>	Specific for Verotoxin II - A subunit (SLT-2a).
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Verotoxins (Shiga-like) antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	Verotoxins antibody (Shiga-like) is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>STORAGE CONDITIONS:</b>	Verotoxins (Shiga-like) antibody can be stored at -20°C, stable for one year.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	Mouse IgG2
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	Anti-Verotoxins (Shiga-like)
-------------------------	------------------------------

### Background

FOR RESEARCH USE ONLY

December 12, 2016