



prosci-inc.com



**ProSci Incorporated**  
12170 Flint Place  
Poway, CA 92064



**Toll Free:** +1 (888) 513 9525  
**Local:** +1 (858) 513 2638  
**Fax:** +1 (858) 513 2692

techsupport@prosci-inc.com

## Bacillus spores Antibody (biotin)

CATALOG NUMBER: 35-112

### Specifications

#### SPECIES REACTIVITY:

#### TESTED APPLICATIONS:

**USER NOTE:** Optimal dilutions for each application to be determined by the researcher.

**IMMUNOGEN:** Purified spores of *B. cereus* + *B. subtilis*

**HOST SPECIES:** Rabbit

### Properties

**PURIFICATION:** Purified IgG fraction of antiserum covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution.

**PHYSICAL STATE:** Liquid

**BUFFER:** Bacillus spores antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

**STORAGE CONDITIONS:** Bacillus spores antibody can be stored at -20°C, stable for one year.

**CLONALITY:** Polyclonal

**CONJUGATE:** Biotin

### Additional Info

**ALTERNATE NAMES:** Anti-Bacillus spores Biotin conjugate

### Background

FOR RESEARCH USE ONLY

December 12, 2016