



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064

HIGH PERFORMANCE ANTIBODIES ... AND MORE

Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Bordetella pertussis Antibody [5113] (FITC)

CATALOG NUMBER: 35-592

Specifications

SPECIES REACTIVITY:	Bacteria
TESTED APPLICATIONS:	ELISA, IF, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Specific for LOS-A of B. pertussis & B. bronchiseptica.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Bordetella pertussis antibody is conjugated with high purity isomer I of fluorescein isothiocyanate. Care is taken to ensure complete removal of any free fluorescein from the final product.
PHYSICAL STATE:	Liquid
BUFFER:	Bordetella pertussis antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
STORAGE CONDITIONS:	Bordetella pertussis antibody can be stored at -20°C, stable for one year.
CLONALITY:	Monoclonal
ISOTYPE:	Mouse IgG3
CONJUGATE:	FITC

Additional Info

ALTERNATE NAMES:	Anti-Bordetella pertussis
-------------------------	---------------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016