



prosci-inc.com



ProSci Incorporated
12170 Flint Place
Poway, CA 92064

HIGH PERFORMANCE ANTIBODIES ... AND MORE

Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

Clostridium difficile toxin Antibody [5155]

CATALOG NUMBER: 35-598

Specifications

SPECIES REACTIVITY:	Bacteria
TESTED APPLICATIONS:	ELISA
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Specific for toxin/toxoid A of Clostridium difficile. Negative vs. other toxins.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Clostridium difficile toxin antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
PHYSICAL STATE:	Liquid
BUFFER:	Clostridium difficile toxin antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
STORAGE CONDITIONS:	Clostridium difficile toxin antibody can be stored at -20°C, stable for one year.
CLONALITY:	Monoclonal
ISOTYPE:	Mouse IgG1
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	Anti-Clostridium difficile toxin
-------------------------	----------------------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016