



prosci-inc.com



**ProSci Incorporated**  
12170 Flint Place  
Poway, CA 92064



**Toll Free:** +1 (888) 513 9525  
**Local:** +1 (858) 513 2638  
**Fax:** +1 (858) 513 2692

techsupport@prosci-inc.com

## Clostridium difficile toxin Antibody [5156]

CATALOG NUMBER: 35-599

### Specifications

<b>SPECIES REACTIVITY:</b>	Bacteria
<b>TESTED APPLICATIONS:</b>	ELISA
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>SPECIFICITY:</b>	Specific for toxin/toxoid B of Clostridium difficile. Negative vs. other toxins.
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Clostridium difficile toxin antibody is purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	Clostridium difficile toxin antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>STORAGE CONDITIONS:</b>	Clostridium difficile toxin antibody can be stored at -20°C, stable for one year.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	Mouse IgG1
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	Anti-Clostridium difficile toxin
-------------------------	----------------------------------

### Background

FOR RESEARCH USE ONLY

December 12, 2016