



## Rb (phospho Ser821) Antibody

CATALOG NUMBER: 47-108

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	WB
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>POSITIVE CONTROL:</b>	1) Cat. No. 1205 - Jurkat Cell Lysate
<b>IMMUNOGEN:</b>	Phosphopeptide region of human Rb containing serine 821.
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	1 mg/ml BSA, 0.05% sodium azide.
<b>STORAGE CONDITIONS:</b>	Rb antibody can be stored at 4°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	p95, Nibrin, NBS, NBS1, P95, Cell cycle regulatory protein p95, ATV, AT-V1, AT-V2
<b>ACCESSION NO.:</b>	NP_002476.2
<b>PROTEIN GI NO.:</b>	33356172
<b>OFFICIAL SYMBOL:</b>	NBN
<b>GENE ID:</b>	4683

### Background

FOR RESEARCH USE ONLY

December 13, 2016