



Rb (phospho Ser811) Antibody

CATALOG NUMBER: 47-107

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
POSITIVE CONTROL:	1) Cat. No. 1205 - Jurkat Cell Lysate
IMMUNOGEN:	Phosphopeptide region of human Rb containing serine 811.
HOST SPECIES:	Rabbit

Properties

PHYSICAL STATE:	Liquid
BUFFER:	1 mg/ml BSA, 0.05% sodium azide.
STORAGE CONDITIONS:	Rb antibody can be stored at 4°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	REPA2, RPA32, RPA34, Replication protein A 32 kDa subunit, Replication factor A protein 2, RP-A p32
ACCESSION NO.:	NP_002937.1
PROTEIN GI NO.:	4506585
OFFICIAL SYMBOL:	RPA2
GENE ID:	6118

Background

FOR RESEARCH USE ONLY

December 13, 2016