



## Giardia lamblia Antibody

CATALOG NUMBER: 35-189

### Specifications

<b>SPECIES REACTIVITY:</b>	Protozoan
<b>TESTED APPLICATIONS:</b>	ELISA, IF, WB
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>SPECIFICITY:</b>	Specific to Intact cysts
<b>IMMUNOGEN:</b>	purified cysts
<b>HOST SPECIES:</b>	Goat

### Properties

<b>PURIFICATION:</b>	Purified IgG fraction of antiserum to a purity of >90%.
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	Giardia lamblia antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
<b>STORAGE CONDITIONS:</b>	Giardia lamblia antibody can be stored at -20°C, stable for one year.
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	Anti-Giardia lamblia
-------------------------	----------------------

### Background

**FOR RESEARCH USE ONLY**

December 12, 2016