



Giardia lamblia Antibody (HRP)

CATALOG NUMBER: 35-191

Specifications

SPECIES REACTIVITY:	Protozoan
TESTED APPLICATIONS:	ELISA, IF, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	Specific to Intact cysts
IMMUNOGEN:	purified cysts
HOST SPECIES:	Goat

Properties

PURIFICATION:	Purified IgG fraction of antiserum covalently coupled with highly purified preparation of horseradish peroxidase (RZ>3). Care is taken to insure adequate conjugation while preserving the maximum enzyme activity. Free enzyme is absent.
PHYSICAL STATE:	Liquid
BUFFER:	Giardia lamblia antibody is in a phosphate saline buffer with bovine serum albumin (10mg/ml) and preserved with 0.002% thimerosal.
STORAGE CONDITIONS:	Giardia lamblia antibody can be stored at -20°C, stable for one year.
CLONALITY:	Polyclonal
CONJUGATE:	HRP

Additional Info

ALTERNATE NAMES:	Anti-Giardia lamblia Peroxidase conjugate
-------------------------	-------------------------------------------

Background

FOR RESEARCH USE ONLY

December 12, 2016