




HMx2 Antibody [1B10]

CATALOG NUMBER: 51-717



Immunohistochemistry of human brain
cortex tissue stained using HMx2
Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC, WB
APPLICATIONS:	HMx2 antibody can be used in Western Blot starting at 0.5 - 2 ug/mL, and immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	HMx2 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG2a,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	HMx2, Homeobox protein HMx2, H6 family homeobox 2, Homeo box (H6 family) 2, Nkx5-2, H6L
ACCESSION NO.:	A2RU54
PROTEIN GI NO.:	152032526
OFFICIAL SYMBOL:	HMx2
GENE ID:	3167

Background

