




CD45RA Antibody [MB1]

CATALOG NUMBER: 70-005



Immunohistochemistry of human spleen tissue stained using CD45RA Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	IHC
APPLICATIONS:	CD45RA antibody can be used in ELISA, Western Blot, and immunohistochemistry starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Tissue Culture Supernatant
PHYSICAL STATE:	Liquid
BUFFER:	0.1% sodium azide
STORAGE CONDITIONS:	CD45RA antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	CD45, Receptor-type tyrosine-protein phosphatase C, Leukocyte common antigen, L-CA, LCA, LY5, B220, T200, CD45R, GP180
ACCESSION NO.:	P08575
PROTEIN GI NO.:	33112650
OFFICIAL SYMBOL:	PTPRC
GENE ID:	5788

Background

FOR RESEARCH USE ONLY

