




IL-6 Antibody, Monoclonal

CATALOG NUMBER: 70-020



Immunohistochemistry of rat spleen tissue
stained using IL-6 Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Rat
TESTED APPLICATIONS:	FACS, IHC, WB
APPLICATIONS:	IL-6 antibody can be used in immunohistochemistry starting at 10 ug/mL.
USER NOTE:	This product ships on dry ice.
IMMUNOGEN:	Purified recombinant rat IL-6.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein G Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2. No stabilizers or preservatives added.
STORAGE CONDITIONS:	IL-6 antibody should be stored at -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG2b
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	IL6, B-cell stimulatory factor 2, BSF-2, BSF2, CTL differentiation factor, Hybridoma growth factor, IFN-beta-2, Interleukin-6, Interferon beta-2, Interleukin BSF-2, IFNB2, B-cell differentiation factor, CDF, HSF, IL-6
ACCESSION NO.:	P05231
PROTEIN GI NO.:	124347
OFFICIAL SYMBOL:	IL6
GENE ID:	3569

Background

FOR RESEARCH USE ONLY

