



PTPN12 Antibody [4G6]

CATALOG NUMBER: 70-151

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IHC, WB
APPLICATIONS:	PTPN12 Monoclonal Antibody can be used for immunohistochemistry (tested on parafin embedded tissues) starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2
STORAGE CONDITIONS:	PTPN12 Monoclonal antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG2a,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	PTPN12, Phosphatase pest, PTPG1, Ptp-g1, PTP-PEST
ACCESSION NO.:	Q05209
PROTEIN GI NO.:	317373522
OFFICIAL SYMBOL:	PTPN12
GENE ID:	5782

Background

FOR RESEARCH USE ONLY

December 13, 2016