



## C1QBP Antibody, Monoclonal

CATALOG NUMBER: 70-253

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	ELISA, FACS, IF, IHC, WB
<b>APPLICATIONS:</b>	C1QBP Monoclonal Antibody can be used for immunohistochemistry (tested on parafin embedded tissues) starting at 1:50 and Western Blot starting at 1:2000 - 1:10000.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>SPECIFICITY:</b>	Recognizes the mature form of gamma C1qR, Mr 33kD and the truncated form lacking aa74-95. Species cross-reactivity: Human.
<b>IMMUNOGEN:</b>	Human monocytes
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Ascites
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	30% glycerol.
<b>STORAGE CONDITIONS:</b>	C1QBP Monoclonal antibody may be stored at 4°C; Avoid freeze-thaw cycles.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	IgG1
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	C1QBP, GC1Q-R, GC1q-R protein, GC1QBP, HABP1, GC1qR, p33, p32, Glycoprotein gC1qBP, Hyaluronan-binding protein 1, SF2p32
<b>ACCESSION NO.:</b>	Q07021
<b>PROTEIN GI NO.:</b>	730772
<b>OFFICIAL SYMBOL:</b>	C1QBP
<b>GENE ID:</b>	708

### Background

FOR RESEARCH USE ONLY

December 13, 2016