



prosci-inc.com



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064



Toll Free: +1 (888) 513 9525  
Local: +1 (858) 513 2638  
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com

## OAT Antibody [AT23A2]

CATALOG NUMBER: 70-287

### Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	IHC
APPLICATIONS:	OAT antibody can be used for staining paraffin embedded tissues starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
HOST SPECIES:	Mouse

### Properties

PURIFICATION:	Protein G Column
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

### Additional Info

ALTERNATE NAMES:	OAT, GACR, HOGA, OKT, Gyrate atrophy, OATASE, Ornithine aminotransferase
ACCESSION NO.:	P04181
PROTEIN GI NO.:	129018
OFFICIAL SYMBOL:	OAT
GENE ID:	4942

### Background

FOR RESEARCH USE ONLY

December 13, 2016