



S100A9 Antibody

CATALOG NUMBER: 70-532

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	IHC, WB
APPLICATIONS:	S100A9 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 10 ug/mL, and Western Blotting starting at 1:200 - 1:2000.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Synthetic peptide
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Antiserum
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	S100A9, 60B8AG, Calgranulin B, CFAG, CAGB, Calgranulin-B, Calprotectin L1H subunit, CGLB, L1Ag, MRP14, NIF, p14, Protein S100-A9, LIAG, MAC387, MRP-14
ACCESSION NO.:	P06702
PROTEIN GI NO.:	115444
OFFICIAL SYMBOL:	S100A9
GENE ID:	6280

Background

FOR RESEARCH USE ONLY

December 13, 2016