




STAT4 Antibody [1C2-1C12]

CATALOG NUMBER: 51-808



Immunohistochemistry of human testis
tissue stained using STAT4 Monoclonal
Antibody.

Specifications

| | |
|-----------------------------|----------------------------------------------------------------------------------------|
| SPECIES REACTIVITY: | Human |
| TESTED APPLICATIONS: | ELISA, IF, IHC, IP, WB |
| APPLICATIONS: | STAT4 antibody can be used in Western Blot, and immunohistochemistry starting at 1:50. |
| USER NOTE: | Optimal dilutions for each application to be determined by the researcher. |
| IMMUNOGEN: | Recombinant protein |
| HOST SPECIES: | Mouse |

Properties

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| PURIFICATION: | Protein A Column |
| PHYSICAL STATE: | Liquid |
| BUFFER: | PBS, pH 7.2 |
| STORAGE CONDITIONS: | STAT4 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles. |
| CLONALITY: | Monoclonal |
| ISOTYPE: | IgG1,k |
| CONJUGATE: | Unconjugated |

Additional Info

| | |
|-------------------------|---------------|
| ALTERNATE NAMES: | STAT4, SLEB11 |
| ACCESSION NO.: | Q14765 |
| PROTEIN GI NO.: | 6226158 |
| OFFICIAL SYMBOL: | STAT4 |
| GENE ID: | 6775 |

Background

