




## STXBP1 Antibody [6D1]

CATALOG NUMBER: 51-811



Immunohistochemistry of human brain,  
Neuropil stained using STXBP1  
Monoclonal Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human, Mouse
<b>TESTED APPLICATIONS:</b>	ELISA, IHC, WB
<b>APPLICATIONS:</b>	STXBP1 antibody can be used in Western Blot, immunohistochemistry starting at 10 ug/mL, immunoprecipitation, and flow cytometry.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	Recombinant protein
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Protein A Column
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2
<b>STORAGE CONDITIONS:</b>	STXBP1 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	IgG2a,k
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	STXBP1, EIEE4, HUNC18, MUNC18-1, N-Sec1, NSEC1, Protein unc-18 homolog 1, Protein unc-18 homolog A, UNC18, Unc18-1, p67, RBSEC1, Unc-18A, Munc18, Neuronal SEC1, Syntaxin binding protein 1, Syntaxin-binding protein 1, UNC18A
<b>ACCESSION NO.:</b>	P61764
<b>PROTEIN GI NO.:</b>	50403646
<b>OFFICIAL SYMBOL:</b>	STXBP1
<b>GENE ID:</b>	6812

