



PTK2 / FAK Antibody

CATALOG NUMBER: 70-377

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat
TESTED APPLICATIONS:	IHC, IP, WB
APPLICATIONS:	PTK2 / FAK Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 1:200, Immunoprecipitation and Western Blotting starting at 1:5000.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant human protein purified from E. coli (His-Fak1. a.a.706-1052).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Affinity Purified
BUFFER:	HEPES with 0.15M NaCl, 0.01% BSA, 0.03% sodium azide, and 50% glycerol
STORAGE CONDITIONS:	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG1
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	PTK2, FADK, Fack1, Focal adhesion kinase, FRNK, FADK 1, Focal adhesion kinase 1, FAK1, p125FAK, PPP1R71, Protein-tyrosine kinase 2, PTK2 protein tyrosine kinase 2, FAK, pp125 FA kinase, Pp125FAK
ACCESSION NO.:	Q05397
PROTEIN GI NO.:	3183518
OFFICIAL SYMBOL:	PTK2
GENE ID:	5747

Background

FOR RESEARCH USE ONLY

December 13, 2016