




CNP (2,3-cyclic nucleotide-3-phosphodiesterase) Antibody

CATALOG NUMBER: 50-105



Western blot of rat brain lysate showing the specific immunolabeling of the ~46k CNP protein.

Specifications

SPECIES REACTIVITY:	Goat, Human, Mouse, Rat, Sheep
TESTED APPLICATIONS:	IHC, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	46
IMMUNOGEN:	Endogenous rabbit 2,3-cyclic nucleotide-3-phospho-diesterase
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Neat serum
PHYSICAL STATE:	Liquid
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	CNPF, CNPI, Cnp1, CNPII, CNP
ACCESSION NO.:	P13233
PROTEIN GI NO.:	51338709
OFFICIAL SYMBOL:	Cnp
GENE ID:	25275

Background

FOR RESEARCH USE ONLY