




Collagen 1, alpha 1 telopeptide Antibody

CATALOG NUMBER: 50-106

Anti-Collagen I α 1, telopeptide



Western blot of rat lung lysate showing specific immunolabeling of the ~140k collagen protein.

Western blot of rat lung lysate showing specific immunolabeling of the ~140 collagen protein.

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat, Sheep
TESTED APPLICATIONS:	IHC, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	140
IMMUNOGEN:	Synthetic peptide from the C-terminal region of human collagen I α 1 telopeptide sequence.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Affinity Purified
PHYSICAL STATE:	Liquid
BUFFER:	100 uL in PBS
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	OI4, Alpha-1 type I collagen,
ACCESSION NO.:	P02452
PROTEIN GI NO.:	296439504
OFFICIAL SYMBOL:	COL1A1
GENE ID:	1277

Background

FOR RESEARCH USE ONLY