




## Nuclear pore Complex Antibody [39C7]

CATALOG NUMBER: 50-132



Immunostaining of HeLa cells with anti-nuclear pore complex antibody (green), and chicken anti-vimentin (red) (Cat. No. 50-264).

### Specifications

<b>SPECIES REACTIVITY:</b>	Human, Mouse, Rat, Yeast
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	Yeast nuclear preparations.
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Total IgG fraction
<b>PHYSICAL STATE:</b>	Liquid
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	IgG1
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	ADIR2, NUP96, NUP196, ADAR2, 98 kDa nucleoporin, Nup98
<b>ACCESSION NO.:</b>	P52948
<b>PROTEIN GI NO.:</b>	308153660
<b>OFFICIAL SYMBOL:</b>	NUP98
<b>GENE ID:</b>	4928

### Background

FOR RESEARCH USE ONLY