




GluR1 Antibody [RH95]

CATALOG NUMBER: 50-119



Western blot of rat hippocampal lysate showing specific immunolabeling of the ~105k GluR1 protein.

Specifications

SPECIES REACTIVITY:	Mouse, Rat
TESTED APPLICATIONS:	WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	100
IMMUNOGEN:	Peptide corresponding to amino acid residues from the N-terminal region of rat GluR1.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein G purified
PHYSICAL STATE:	Liquid
BUFFER:	100 uL in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 ug per mL BSA and 50% glycerol.
CLONALITY:	Monoclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	GluA1, gluR-A, Glur1, AMPA-selective glutamate receptor 1, GluR-1
ACCESSION NO.:	P19490
PROTEIN GI NO.:	97536283
OFFICIAL SYMBOL:	Gria1
GENE ID:	50592

Background

FOR RESEARCH USE ONLY

