




Mitochondria Antibody [113-1]

CATALOG NUMBER: 50-124



Western blot of HeLa lysate showing specific immunolabeling of the ~60k mitochondrial protein.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	IF, IHC, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	60
IMMUNOGEN:	Human cell homogenate
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Total IgG fraction
PHYSICAL STATE:	Liquid
CLONALITY:	Monoclonal
ISOTYPE:	IgG1
CONJUGATE:	Unconjugated

Additional Info

Background

FOR RESEARCH USE ONLY