




## Testicular Receptor 4 (TR4) Antibody

CATALOG NUMBER: 50-147



Western blot of rat cerebellum lysate showing specific immunolabeling of the ~64k TR4 protein.

Western blot of rat cerebellum lysate showing specific immunolabeling of the ~64k TR4 protein.

### Specifications

<b>SPECIES REACTIVITY:</b>	Mouse, Rat
<b>TESTED APPLICATIONS:</b>	
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>PREDICTED MOLECULAR WEIGHT:</b>	64
<b>IMMUNOGEN:</b>	Fusion protein of N-terminal peptide from mouse TR4.
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Neat Serum
<b>PHYSICAL STATE:</b>	Liquid
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	Tr4, TAK1, mKIAA4145, Mtr2r1, Tak1, Tr4, Orphan nuclear receptor TAK1,
<b>ACCESSION NO.:</b>	P49117
<b>PROTEIN GI NO.:</b>	1351191
<b>OFFICIAL SYMBOL:</b>	Nr2c2
<b>GENE ID:</b>	22026

### Background

FOR RESEARCH USE ONLY