




## Glial Fibrillary Acidic Proetin (GFAP) Antibody

CATALOG NUMBER: 50-109



Western blot of rat cortex lysate showing specific immunolabeling of the ~50k GFAP protein.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>PREDICTED MOLECULAR WEIGHT:</b>	50
<b>IMMUNOGEN:</b>	Recombinant and purified bovine GFAP.
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Neat Serum
<b>PHYSICAL STATE:</b>	Liquid
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	GFAP
<b>ACCESSION NO.:</b>	Q28115
<b>PROTEIN GI NO.:</b>	143811396
<b>OFFICIAL SYMBOL:</b>	GFAP
<b>GENE ID:</b>	281189

### Background

FOR RESEARCH USE ONLY