




Synaptobrevin (VAMP) Antibody [SP10]

CATALOG NUMBER: 50-142



Western blot of rat brain lysate showing specific immunolabeling of the ~16k synaptobrevin protein.

Specifications

SPECIES REACTIVITY:	Bovine, Hamster, Human, Mouse, Rat
TESTED APPLICATIONS:	IHC, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	16
IMMUNOGEN:	Synaptic immunoprecipitate (crude) from human brain.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	IgM Purified
PHYSICAL STATE:	Liquid
BUFFER:	100 uL in 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 ug per mL BSA and 50% glycerol.
CLONALITY:	Monoclonal
ISOTYPE:	igM
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	SYB1, VAMP-1, SYB1, Synaptobrevin-1
ACCESSION NO.:	P23763
PROTEIN GI NO.:	135093
OFFICIAL SYMBOL:	VAMP1
GENE ID:	6843

Background

FOR RESEARCH USE ONLY

