




## Periostin pan Antibody

CATALOG NUMBER: 50-139



Western blot of rat lung lysate showing specific immunolabeling of the ~93 kDa periostin protein triplet.

### Specifications

<b>SPECIES REACTIVITY:</b>	Chicken, Human, Mouse, Rat
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>PREDICTED MOLECULAR WEIGHT:</b>	93
<b>IMMUNOGEN:</b>	peptide from fasciclin domain 1 of mouse periostin.
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Affinity Purified
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	100 uL in PBS
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	PN, PLF, Osf2, peri, OSF-2, AI747096, A630052E07Rik, Osteoblast-specific factor 2, PN
<b>ACCESSION NO.:</b>	Q62009
<b>PROTEIN GI NO.:</b>	46576895
<b>OFFICIAL SYMBOL:</b>	Postn
<b>GENE ID:</b>	50706

### Background

FOR RESEARCH USE ONLY