




TDP43 Antibody

CATALOG NUMBER: 50-145



Western blot of neonatal rat brain lysate showing specific immunolabeling of the ~43 kDa TDP43 protein.

Specifications

SPECIES REACTIVITY:	Human, Rat
TESTED APPLICATIONS:	IHC, WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	43
IMMUNOGEN:	Synthetic peptide C-terminal to the caspase-cleavage site (between D219 and V220) of human TDP-43.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Affinity Purified
PHYSICAL STATE:	Liquid
BUFFER:	100 uL in PBS + 30% glycerol
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	ALS10, TDP-43, TDP43
ACCESSION NO.:	Q13148
PROTEIN GI NO.:	20140568
OFFICIAL SYMBOL:	TARDBP
GENE ID:	23435

Background

FOR RESEARCH USE ONLY