




Testicular Receptor 2 (TR2) Antibody

CATALOG NUMBER: 50-146



Western blot of rat testes lysate showing specific immunolabeling of the ~64k TR2 protein.

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat
TESTED APPLICATIONS:	WB
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
PREDICTED MOLECULAR WEIGHT:	64
IMMUNOGEN:	Fusion protein from N-terminal peptide of mouse TR2.
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Neat Serum
PHYSICAL STATE:	Liquid
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	TR2, 80.3, Eenr, Tr2-11, 4831444H07Rik, Tr2, Orphan nuclear receptor TR2, mTR2
ACCESSION NO.:	Q505F1
PROTEIN GI NO.:	341942228
OFFICIAL SYMBOL:	Nr2c1
GENE ID:	22025

Background

FOR RESEARCH USE ONLY