




Tuj1 Antibody

CATALOG NUMBER: 48-004



Immunohistochemistry of human testis tissue stained using Tuj1 Antibody.

Specifications

SPECIES REACTIVITY:	Human, Mouse, Rat
TESTED APPLICATIONS:	ICC, IHC, WB
APPLICATIONS:	Tuj1 antibody can be used in Western Blot, immunohistochemistry starting at 1:50, immunocytochemistry, and flow cytometry.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Peptide corresponded to a region of the Tuj 1 gene product shared between the rat (AAM28438, NCBI) and human (AAL28094, NCBI) gene products.
HOST SPECIES:	Chicken

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.2, 0.02% sodium azide
STORAGE CONDITIONS:	Tuj1 antibody may be stored at 4 °C for short-term only. As with all antibodies avoid freeze/thaw cycles. Store at -20 °C for long term. Aliquots are stable for at least 12 months
CLONALITY:	Polyclonal
ISOTYPE:	IgY
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	TUBB3, CDCBM, Beta III Tubulin, Class III beta-tubulin, TUBB4, Tubulin, beta 3, Tubulin beta-III, Tubulin beta-3 chain, Tubulin beta-4 chain, Tubulin, beta 3 class III, CFEOM3A
ACCESSION NO.:	Q13509
PROTEIN GI NO.:	20455526
OFFICIAL SYMBOL:	TUBB3
GENE ID:	10381

FOR RESEARCH USE ONLY

December 13, 2016