




## BAIAP3 Antibody

CATALOG NUMBER: 48-651



Immunohistochemistry staining of BAIAP3  
in brain cortex tissue using BAIAP3  
Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	ELISA, IHC
<b>APPLICATIONS:</b>	BAIAP3 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 1:500, immunofluorescence, and immunoprecipitation.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	BAIAP3 antibody was raised against amino acids 1171 to 1187 of BAIAP3 (Human).
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Protein G Column
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, 0.09% sodium azide.
<b>STORAGE CONDITIONS:</b>	BAIAP3 antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Polyclonal
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	BAIAP3, BAI-associated protein 3, BAI1-associated protein 3, BAP3, KIAA0734
<b>ACCESSION NO.:</b>	O94812
<b>PROTEIN GI NO.:</b>	51701298
<b>OFFICIAL SYMBOL:</b>	BAIAP3
<b>GENE ID:</b>	8938

### Background

FOR RESEARCH USE ONLY

