




## TTRAP Antibody

CATALOG NUMBER: 48-775



Immunohistochemistry staining of TTRAP  
in kidney tissue using TTRAP Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Gibbon, Gorilla, Human
<b>TESTED APPLICATIONS:</b>	IHC, WB
<b>APPLICATIONS:</b>	TTRAP antibody can be used in immunohistochemistry starting at 1:50, and flow cytometry starting at 1:1 - 1:10.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	TTRAP antibody was raised against amino acids 11 - 28 of TTRAP (Human).
<b>HOST SPECIES:</b>	Rabbit

### Properties

<b>PURIFICATION:</b>	Immunoaffinity Chromatography
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, 0.02% sodium azide.
<b>STORAGE CONDITIONS:</b>	Store TTRAP antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Polyclonal
<b>ISOTYPE:</b>	IgG
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	TDP2, 5-Tyr-DNA phosphodiesterase, AD022, EAP2, EAPII, ETS1-associated protein 2, ETS1-associated protein II, HTDP2, RP1-30M3.3, TTRAP, DJ30M3.3, Tyr-DNA phosphodiesterase 2
<b>ACCESSION NO.:</b>	O95551
<b>PROTEIN GI NO.:</b>	67462008
<b>OFFICIAL SYMBOL:</b>	TDP2
<b>GENE ID:</b>	51567

### Background

FOR RESEARCH USE ONLY

