




BRAF Antibody

CATALOG NUMBER: 48-778



Immunohistochemistry staining of BRAF in testis tissue using BRAF Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	IF, IHC
APPLICATIONS:	BRAF antibody can be used in immunohistochemistry starting at 1:50.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	pThr598/Ser601
IMMUNOGEN:	BRAF antibody was raised against a modified peptide of BRAF (Human).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	Phosphate-buffered solution, pH 7.2, 0.09% sodium azide, 50% glycerol.
STORAGE CONDITIONS:	Store BRAF antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
ISOTYPE:	IgG
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	BRAF, 94 kDa B-raf protein, B-raf, BRAF1, B Raf, Proto-oncogene B-Raf, p94, B-RAF1, NS7, ONCOGENE BRAF, RAFB1
ACCESSION NO.:	P15056
PROTEIN GI NO.:	50403720
OFFICIAL SYMBOL:	BRAF
GENE ID:	673

Background

