




NTSR2 Antibody

CATALOG NUMBER: 51-245



Immunohistochemistry staining of NTSR2
in adrenal tissue using NTSR2
Antibody.

Specifications

SPECIES REACTIVITY:	Gibbon, Human
TESTED APPLICATIONS:	ELISA, IHC
APPLICATIONS:	NTSR2 antibody can be used in immunohistochemistry starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	BLAST analysis of the peptide immunogen showed no homology with other human proteins.
IMMUNOGEN:	NTSR2 antibody was raised against amino acids from the 3rd extracellular domain of NTSR2 (Human).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.1% sodium azide.
STORAGE CONDITIONS:	NTSR2 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	NTSR2, Neurotensin NT2 receptor, Neurotensin receptor, type 2, NTR2, NT2 receptor, Neurotensin receptor 2, Neurotensin receptor type 2, NT-R-2
ACCESSION NO.:	O95665
PROTEIN GI NO.:	311033404
OFFICIAL SYMBOL:	NTSR2
GENE ID:	23620

Background

