




Integrin Alpha 5 Antibody [10F6]

CATALOG NUMBER: 51-385



Immunohistochemistry of human heart
capillaries stained using Integrin Alpha 5
Monoclonal Antibody.

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, FACS, IHC, WB
APPLICATIONS:	Integrin Alpha 5 antibody can be used in Western Blot, and immunohistochemistry starting at 1 - 2 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	Recombinant protein
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Ascites
PHYSICAL STATE:	Liquid
BUFFER:	0.03% sodium azide
STORAGE CONDITIONS:	Integrin Alpha 5 antibody should be stored long term (months) at -20 °C and short term (weeks) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG2a
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	ITGA5, FNRA, Integrin alpha-5, Integrin Alpha 5, Integrin alpha-F, VLA-5, CD49e, CD49e antigen, VLA5A
ACCESSION NO.:	P08648
PROTEIN GI NO.:	23831237
OFFICIAL SYMBOL:	ITGA5
GENE ID:	3678

Background

