




SH3BP5 Antibody [2B3]

CATALOG NUMBER: 51-543



Immunoperoxidase of monoclonal antibody to SH3BP5 on formalin-fixed paraffin-embedded human colon adenocarcinoma (antibody concentration 1 ug/mL).

Specifications

SPECIES REACTIVITY:	Human
TESTED APPLICATIONS:	ELISA, IF, IHC, WB
APPLICATIONS:	SH3BP5 antibody can be used in Western Blot, and immunohistochemistry starting at 5 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
IMMUNOGEN:	SH3BP5 (NP_004835, 4 a.a. ~ 113 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
HOST SPECIES:	Mouse

Properties

PURIFICATION:	Protein A Column
PHYSICAL STATE:	Liquid
BUFFER:	PBS, pH 7.4.
STORAGE CONDITIONS:	Store SH3BP5 antibody at -20 °C. Aliquot to avoid freeze/thaw cycles.
CLONALITY:	Monoclonal
ISOTYPE:	IgG1,k
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	SH3BP5, SH3 domain-binding protein 5, SH3BP-5, SAB
ACCESSION NO.:	O60239
PROTEIN GI NO.:	124028536
OFFICIAL SYMBOL:	SH3BP5
GENE ID:	9467

Background

