




## G3BP Antibody [2F3]

CATALOG NUMBER: 51-549



Immunohistochemistry of human breast tissue stained using G3BP Monoclonal Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	ELISA, IF, IHC, RNAi, WB
<b>APPLICATIONS:</b>	G3BP antibody can be used in Western Blot, and immunohistochemistry starting at 1:200 - 1:400.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	G3BP (AAH06997, 214 a.a. ~ 303 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Protein A Column
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2
<b>STORAGE CONDITIONS:</b>	Store G3BP antibody at -20 °C. Aliquot to avoid freeze/thaw cycles.
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	IgG1,k
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	G3BP1, HDH-VIII, G3BP-1, G3BP, GAP binding protein, HDH VIII
<b>ACCESSION NO.:</b>	Q13283
<b>PROTEIN GI NO.:</b>	14916572
<b>OFFICIAL SYMBOL:</b>	G3BP1
<b>GENE ID:</b>	10146

### Background

